



SKR-AB00C

All-in-One Computer Cart



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

TECHNICAL PARAMETERS

Size: 850*520*920mm

QUALITY CONTROL

COMPUTER HOLDER

Suitable for almost all types of PC



DOUBLE SIDED TELESCOPIC HOLDER

Double sided telescopic holder for mouse, left or right hand operating is both available



WHEELS

Casters is using medical silent material, double sided casters, stable and noiseless. Front 2 casters with brake, lockable at 360 degree, rear 2 casters free and swivel to any direction.



BASE

Base is ABS covered high quality metal structure, featured with both stability and easy-maintenance.



ABS MATEIRAL TOP

Easy to clean and endurable to scratch. Cover is free to open, laptop can be put in easily



MATERIAL

Metal structure material comes from Baosteel from Fortune Global 500





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-1802123190 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT: 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.