


GOT YOU

BRAND IS OUR ETERNAL PURSUIT

SKB038-2

MEDICAL CARE SPECIALIST



Our pricing, no middlemen
make a difference



Factory



Network



Customer



0-75°

Back-rest

Allows the resident to get out of bed with confidence



620-860_{mm}

Height Adjustment

All four side rails are raised entirely to provide patient security



Detail display



Gas spring
Manual operate the
back rest position



Crank
Adjust trolley height



5" castors with brake
Deluxe mute wheel

The shape is similar and the quality is different

The peers have been imitating all the time, but we have never been surpassed.

Similar in shape and different in process

1 Line spray exquisite workmanship High hardness, high wear resistance



✘ Manual spray, rough technology, easy to fall off the surface.

Similar in shape and different in process

2 Big ABS rail Sturdy, durable, soft treatment



✘ Cheap accessories, easy to break, and have a safety hazard.

Introduction

Size: 1960*830*620-860mm

Gas spring adjust back rest 0-75°

With mattress

5" castors with brake


Powder coated steel structure

Crank adjust trolley height

One pair of big ABS rail

Multiple program selection, full automatic intelligent care

Reliable quality and good faith service



About us

Medical sickbed expert

Jiangsu Saikang Medical Equipment Co.,Ltd. has been established since 2002. It is located in Zhangjiagang City, along the Economic circle of Yangtse River Delta, enjoying the highly developed economic expansion and convenient transportation access.

FDA



CE

