



SKH098-3

Cabinet



SAIKANG / Medical


LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098

TECHNICAL PARAMETERS

Size: 900 x 500 x 1800mm

QUALITY CONTROL



12_{mm}

THIN FRAME DESIGN
Thin Frame Design(12mm), bigger space, easier to take and storage the medicine, instruments or files.



DESIGN FOR HANDLE
Handle designed inside the door, avoid possible scratch with user



THICKNESS OF MATERIA
0.6mm steel thickness, 0.12mm powder coating, while market thickness is 0.5mm in total.


0.6_{mm} | 0.12_{mm} | 0.5_{mm}



WORLD BRAND MATERIAL SOURCE
Baosteel from Fortune Global 500



LOCK DESIGN
Lock column lower side of door, prevent hurting users.



POWDER COATING
Process Epoxy Painting, ASTM Testing Anti-Baterial, Paint Thickness 0.12mm, Brightness 60°, Paint Can Resist 50kg Impact

0.12_{mm} | 60° | 50_{kg}



KNOCK DOWN PACKAGE
Knock-down packing only need 0.25CBM, 1/3 of assembled packing (0.76CBM), package with foam and other materials to protect fragile parts like glass, etc.

0.25_{CBM}



COLOR OPTIONS
White color, Blue color other color optional

RICH PARTS

LEADING BRAND OF CHINA MEDICAL FURNITURE



Label socket



Adjust Hole



Drawer Rail



Hinge



Latch



Mirror





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.