



A100-4A

Electric
Operating
Table

SKL-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098

TECHNICAL PARAMETERS

External L*W*H	2060*500*(710-1010)mm
Reverse Trendelenburg	0-25°
Trendelenburg	0-25°
Lateral Tilt Left/right	0-15°
Head Plate up/down	0-40°/0-90°
Back Plate up/down	0-80°/0-20°
Leg Plate up/down	0-15°/0-90°
Kidney Bridge up	110mm
Horizontal Sliding C type Arm	350mm

TECHNICAL CONFIGURATION

Bed Frame	1pc
Mattress	1set
Anesthesia Screen Frame	1pc
Shoulder Holder with Cushion	2pcs
Arm Holder with Cushion	2pcs
Leg Support with Cushion	2pcs
Foot Support with Cushion	2pcs
Hand Controller	1pc
Kidney Bridge Handle	1pc

ELECTRONICS

Battery Standby Time(h)	80
Power Supply	220V 50HZ 350VA
Power Control Circuit	12V
Plug	as required

QUALITY CONTROL

MATERIAL

The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned



POWDER SUPPLY

Internal battery maintaining power supply for 50-80 operations and with low battery warning function

50-80



REMOTE PANEL & COLUMN PANEL

Leg Plate With Imported Gas Spring For Easy Control.



BRAKING SYSTEM

Pedal Brake, Easy To Move And Lock.




WHEEL

6" castors import from Germany



HYDRAULIC SYSTEM

Brand hydraulic system (German Hoerbiger) make table with low noise, high safety and reliability



LEG PLATE

Leg Plate With Imported Gas Spring For Easy Control.



PROFESSIONAL FUNCTIONS

LEADING BRAND OF
CHINA MEDICAL FURNITURE

Lowest Position of
Table Top

710mm



Highest Position of
Table Top

1010mm



Back Plate
Upward

0°-80°



Back Plate Downward

0°-20°



Reversed
Trendelenburg

0°-25°



Trendelenburg

0°-25°




Head Plate Upward

0°-40°



Head Plate Downward

0°-90°



Right Lateral Tilt

0-20°



Left Lateral Tilt

0-20°



Horizontal Sliding

350mm



Leg Plate Rest
-90°~15°

RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



Armrest



Anesthesia Frame



Leg Crutches



Shoulder Support



Body Supporter





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

UK Healthward International Limited
Owner(s) Address-DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.