



A98-3

Electric
Obstetric Bed


SKL-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

QUALITY CONTROL




MATERIAL

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean



BED SURFACE

Bed surface has PP cover, made by virgin material. The PP cover is detachable, and easy to clean.




GUARDRAIL

Handrail is rotatable. when you don't use it, you can put down the guardrail.



MATTRESS

The mattress adopts high density memory sponge and is once formed. Mattress is movable.



FOOT SECTION

Foot section can be retracted into the bed, convenient for gynecological and delivery surgical.




CASTER

5inch central controlled caster for moving obstetric table.
TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time

ACCESSORY

Equipped with stainless steel sewage basin and other various accessories



RICH PARTS



Handle



Lithotomy rod



Foot support



Hand controller



Stainless Steel Basin



Caster

DIMENSION PARAMETER

External Length	2050mm
External Width	980mm
Table Height Min.	580mm
Table Height Max.	950mm

TECHNICAL PARAMETERS

Back-rest Adjustment	0-70°
Reverse Trendelenburg	0-5°
Trendelenburg	0-12°
Height Adjustment	580-950mm

STANDARD ACCESSORIES

Pull Handle	2pcs
Leg Support	2pcs
Foot Support	2pcs
Hand Remote Controller	1pc
Stainless Steel Basin	1pc
Instruction Manual	1pc



INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Video

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.