



A99-8

Electric Gynecological
Exam Couch

SKF-C



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

DIMENSION PARAMETER

External Length	1450mm
External Width	600mm
Table Height Min.	580mm
Table Height Max.	1000mm

TECHNICAL PARAMETERS

Height Of Table Top	580mm
Back board upward	70°
Back board downward	20°
Seat board upward	40°
Main Voltage AC	220V±22V,50Hz±1Hz
Patient Weight Capacity	185kg

STANDARD ACCESSORIES

Pull Handle	2pcs
Knee Curtch	2pcs
Pedal Remote Controller	1pc
Stainless Steel Basin	1pc
Instrution Manual	1pc

QUALITY CONTROL




MATERIAL

The base seat and height adjustment system are made of powder coated steel, high gloss, anti corrosion, easy to clean.




PEDAL REMOTE CONTROLLER

Famous brand motor, Various actions are controlled by remote controller.



PULL HANDLE

The hand bar can be dismantled and movable.



STAINLESS STEEL BASIN

Equipped with stainless steel sewage basin and other various accessories.



MATTRESS

The mattress adopts high density memory sponge and is once formed, with no seam.



KNEE CURTCH

Double using leg support, can be used as boatstretcher.

PROFESSIONAL FUNCTIONS

Lowest Position of
Table Top **580mm**



Highest Position of
Table Top **1000mm**



Back Board Upward
70°



Back Board Downward
20°



Seat board upward
40°



RICH PARTS



Pull Handle



Knee Curttch



Pedal Remote Controller



Stainless Steel Basin



Wheel



INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.