



CT2k

PEDIATRIC MANUAL BED

Registration code: SK-C2




SAIKANG / Medical

Healthier
your life

CT2k

Pediatric manual bed



Technical parameters

▪ Size (LxW):	1800x730 mm
▪ Height of use base:	500 mm
▪ Mattress platform:	1580x660 mm
▪ Back-rest adjustment:	0-75° (±10°)
▪ Knee-rest adjustment:	0-35° (±10°)

Technical configuration

<input checked="" type="checkbox"/> Manual crank	2 pcs.
<input checked="" type="checkbox"/> Anti-slip bed feet cover	4 pcs.
<input checked="" type="checkbox"/> Aluminium alloy side rails with extra protection	1 set
<input checked="" type="checkbox"/> Aluminium alloy and wood bed ends with safe lock	1 set
<input checked="" type="checkbox"/> IV pole	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs.
<input checked="" type="checkbox"/> Drainage hook	2 pcs.
<input checked="" type="checkbox"/> Shoes holder	1 set



SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 📞 +86 18021231901 · 🕒 24/7
 ✉ sales@chinasaikang.com
 🌐 www.saikangmedical.com



Company video



General catalogue