



SKH018

ANESTHESIA TROLLEY

Registration code: SKM-B



SKH018

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 600×400×850 mm

Technical configuration

- 4" Covered Castors 4pcs



Material

Made of the high-strength stainless steel, resistance of water, acid, alkali, fire and easy clean, not easily polluted by reagents, blood, drugs etc., satisfy the daily disinfectant and infection control with long life.



Table Top

Guardrail on three side to prevent the falling during moving.




Drawer

Upper drawer with plastic handle prevents cutting, one lock is equipped with two keys.



Cabinet

Lower cabinet with large capacity, plastic handle with lock, easy to storage large volume items.



Castors

4 nylon wheels Ø 100 mm with non-marking tires and plain bearings, two with brakes.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue