



# SKH048

## OXYGEN BOTTLE TROLLEY

Registration code: SKM-B



**SAIKANG** / Medical

Healthier  
your life

# SKH048

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 1435×360×1040 mm

### Technical configuration

- 6" TPP Rubber Rastor 2pcs
- 3" Covered Castors 2pcs



#### Design

U-shaped groove design, the radian meets the requirements of the oxygen bottle for safe transportation.



#### Steels | Welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedding.




#### Dimension

Size: 1435×360×1040mm, Can load 40L oxygen cylinders, the body is light and easy to move.



#### Galvanized Chain

4mm galvanized chain, operational safety.



#### Handle

Pushing handle, easy to move the trolley.



#### Thickening Plate

Cold rolled steel plate, effective compression, bearing weight, not easy to damage, long service life.



#### 6" TPP Rubber Rastor

6" TPP rubber castor, Strong, durable, labor-saving, safty in use.



#### 3" Covered Castors

3" Covered Castors, each with brakes, easy to install, stable and smooth.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue