



SKR-MT625

MEDICINE TROLLEY

Registration code: SKM-B



SAIKANG / Medical

Healthier
your life

SKR-MT625

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 625×475×930 mm

Technical configuration

<input checked="" type="checkbox"/> Dust Baskets	1pc
<input checked="" type="checkbox"/> Utility Container	1pc
<input checked="" type="checkbox"/> Central Lock	1set
<input checked="" type="checkbox"/> Trash Can	2pcs
<input checked="" type="checkbox"/> Castor	4pcs



Powder Coated Column

Make the structure of the trolley more stable and firm.



Drawers

1-4 layer small drawer, fifth floor middle drawer, sixth floor large drawer, can put different sizes of medicine boxes and other medical materials, and equipped with central control lock.



Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley.



Trash Can

There are two ABS dirt buckets of the car body, different colors can distinguish different kinds of dirt.



File Box

Side attached file box, for storage of paper work.



Heavy Duty Construction

Wheels assembled on the same vertical line as trolley four columns, more stable.



Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.



Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



Drawer Card Holder

card insert from downside, easy keep clean.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue