



SKR007

ABS TROLLEY

Registration code: SKM-B



SAIKANG / Medical

Healthier
your life

SKR007

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 730×480×900 mm

Technical configuration

- ☑ 3" covered castors, two with brakes 4 pcs.
- ☑ Plastic bucket (Φ230×200mm) 1 pc.
- ☑ Stainless steel basin (Φ300mm) 1 pc.
- ☑ Φ25, Φ13mm stainless steel round pipe



Structural design

The structure is compact and practical; the use of high-quality thick stainless steel and ABS material, safe and non-toxic, easy to clean, smooth surface, not easy to corrode, and long service life.




Table Top

The countertop is high-quality ABS. The countertop is equipped with three-sided stainless steel fence to prevent items from falling.



Drawer

The drawer is equipped with high-quality silent slide rails, which can be drawn smoothly and smoothly, which effectively solves the problem of easy sliding out or blocking during the drawing process of the drawer, with 3pcs*3 pcs plastic separators.



Application

Stainless steel & ABS utility trolley, Multifunctional nursing trolley, medical Trolley.



Castor

High-quality medical silent universal castors, flexible and light to drive, smooth steering, wear-resistant, with diagonal brake, can be easily braked in any state.



Basin/Bucket

The product adopts a dirt basin/bucket design, the dirt can be put into the basin/bucket, it is convenient to disassemble and assemble, and the basin/bucket can be taken out directly.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue