



# X26

Electric  
Examination Table

SKD-A



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098



Back-rest

Height Adjustment

## TECHNICAL PARAMETERS

External L*W*H	1810*680* (570-860) mm
Back-rest Adjustment	0-75°

## TECHNICAL CONFIGURATION

Pillow	1 pc
Motor	1 pc
Mechanical Lever to Adjust Back Rest	1 pc
2" Casters	2 pcs
Hand Controller	1 pc
Paper Roll Holder	1 pc

## ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
Plug	as required (hospital grade)

## QUALITY CONTROL

### LEATHER

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance

100,000<sub>times</sub>



### SPONGE

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³

4<sub>cm</sub> | 30<sub>kg/m³</sub>



### PAINTING

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

0.12<sub>mm</sub> | 60° | 50<sub>kg</sub>



### METAL MATERIAL

Baosteel from Fortune Global 500  
Panasonic Robotic ensure 360° full smooth welding

500 | 360°



### ELECTRIC SYSTEM


Taiwan motor, Noise under 50dB, with lifetime 20,000times movement

50<sub>dB</sub> | 20,000<sub>times</sub>




### BACK ADJUST

Back adjust can choose ratchet or motor



## RICH PARTS

LEADING BRAND OF  
CHINA MEDICAL FURNITURE



1.2-1.5 mm Thickness powder coated steel frame,  
150 kg bear loading




2 pcs 2' Hidden castors, 2 pcs 2' universal wheels



Hand controller; keep qualified after 2000 times  
use



Engineering ABS material plug, foot pad



Paper roll holder  
Foot bar system  
Motor to adjust the back section





## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

### JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.chinasaikang.com](http://www.chinasaikang.com)

#### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

#### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

#### SAIKANG (UK Branch)

Uk Healthward International Limited  
Owner(s) Address DEPT 114,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

#### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
PO Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

#### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
X22-S03-England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.