



X27


Electric Examination Table

SKD-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098



Back-rest

Waist Section up

Height Adjustment

TECHNICAL PARAMETERS

External L*W*H	1970*690* (530-970) mm
Back-rest Adjustment	0-75°

TECHNICAL CONFIGURATION

Pillow	1pc
Motor	2pcs
Gas Spring	2pcs
4" Casters	4pcs
Hand Controller	1pc
Paper Roll Holder	1pc

ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
Plug	as required (hospital grade)

QUALITY CONTROL

LEATHER

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance

100,000_{times}

SPONGE

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³

4_{cm} | 30_{kg/m³}

PAINTING

1 | Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

0.12_{mm} | 60° | 50_{kg}

METAL MATERIAL

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full smooth welding

500 | 360°

ELECTRIC SYSTEM

Taiwan motor, Noise under 50dB, with lifetime 20,000times movement

50_{dB} | 20,000_{times}

BACK ADJUST


Back adjust can choose ratchet or motor

RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



1.2-1.5 mm Thickness powder coated steel frame,
150 kg bear loading




4 pcs 2" Hidden castors, 4 pcs 2" universal wheels



Hand controller; keep qualified after 2000 times
use



Engineering ABS material plug, foot pad



Paper roll holder
Foot bar system
Motor to adjust the back section





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.