



Z0n

Electric Bath Bed



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

QUALITY CONTROL




HAND CONTROL

Use the special imported waterproof electric motor, making the height adjustment smoothly.

OPTIONAL: Hydraulic pump height adjustment

STEELS | PAINTING

The main frame used high quality S.S. with the Mirror surface polishing



BATH

Made of high toughness and environmentally friendly kayaking materials, the materials pass REACH, ROHS, EN71 and other standards; high density EVA in the pad; soft and comfortable, high temperature/cold cold (+80°-10°), not easy to change, not easy to age; Detachable cleaning and disinfection

WHEEL

Advanced protection universal brake wheel, easy to move, simple to operate, reliable performance



GUARDRAIL

The guardrail is 270mm high, and the guardrails on both sides can be rotated by 180°. The bather can be transferred directly from the bed to the bath; easy to move and move.





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)


Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.