



SK-ET75077A

Emergency Trolley



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

QUALITY CONTROL

MATERIAL

Steel and aluminum structure inside, ABS and antibacterial materials outside, easy to clean.



DRAWER

Drawer insize big & deeper!
Inner with partitions can be well organized



EXTENDABLE WORKTOP

Laterally extendable worktop



TABLE-BOARD

With handles on both sides, surface with transparent soft glass



IV POLE

Movable IV pole can exchange position with defibrillator shelf



DEFIBRILLATOR SHELF

The plate surface is made of virgin ABS material. The position of the IV pole and the defibrillator shelf can be exchanged for easy use.



POWER OUTLET

The power socket is movable and detachable, Easy to use, Located above the CPR board to convenient the use of oxygen cylinders



OXYGEN BOTTLE HOLDER

Inserted oxygen bottle holder, belts fixed on trolley main body



TRASH CAN

Two ABS trash cans with lid



WHEELS

TPR tire no worn out after running 30KM, save anti-winding hard shell, United forming without bolts
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500times



DEBRIS BASKET

With two debris baskets on both sides. The big one is for two sharp containers and the small one is for sanitizer





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-1802123190 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.