




SKB041-6

Patient
Transportation
Trolley



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098



Back-rest

Height Adjustment

MAIN TECHNICAL PARAMETERS

Size	1930*660*(510-850)mm
Back-rest	0-70°(±10°)

TECHNICAL CONFIGURATION

Crank System	1 set
Back Rest Air Spring	1 pc
Plastic Side Rail	1 set
6" Castors , Central Locking	4 pcs
Directional Wheel	1 pc
Directional Pedal	2 pcs
IV Pole Previsions	2 pcs
IV Pole	1 set
Oxygen Cylinder Holder	1 pc
Central Locking Break Pedal	4 pcs
Transfer Mattress	1 pc
Safe Lock for Side Rail	2 pcs

QUALITY CONTROL

GUARDRAIL REVERSIBLE

Side-rail can be Fixed Horizontal Level. Fixable at horizontal level, where patient's arm can be placed on for the administration of IV or other treatment.



T SHAPE DESIGN

T-shaped design allows tubes to pass through side-rail for improved fluid drainage.



DEGREE INDICATE

Scale on guardrail, shows the backrest degree.



ONE DIRECTION MOVING

The side-rail can be laid out one a bed next to the stretcher to act as a transfer board, allows medical personnel transfer patient from stretcher to bed more efficiently and easier.



CENTER WHEEL

Can be engaged to facilitate direction control, enabling users to navigate the stretcher easily in hectic and dynamic environments.

RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



Adjustable vertical oxygen bottle, load capacity:15KG



Stainless steel handle with plastic handle, more durable and beautiful



6 inch controlled castors



Full SS IV pole, bearing 15kg



Water proof mattress with 2 belts





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.