



# SKB1C02-1

Stair Stretcher

SKW-B2



**SAIKANG** / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE  
Stock Code: 870098

## DIMENSION PARAMETER

Size: 131\*50\*155cm

## QUALITY CONTROL



### STRUCTURAL DESIGN

Scientific design, ergonomics, simple structure, simple operation, reliable and secure. Widely used in the home for the old people, school security, hospitals, hotels etc for emergency evacuation.



### FOLDABLE

Foldable stretcher, easy for carrying and storage



### FUNCTION

Backrest height adjustable, meet different users



### WHEELS

With 4 wear-resistant rubber wheels, the chair can use as a wheelchair on the ground.

### MECHANICAL FUNCTION

Can climb up and down stairs smoothly and steadily by rubber track. Can move on the ground as manual wheelchair. With Strong grip track, prevent glide or fall down when stop on the stairs.



### POWDER COATING


11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

0.12<sub>mm</sub> | 60° | 50<sub>kg</sub>



# RICH PARTS

LEADING BRAND OF CHINA MEDICAL FURNITURE




Folding Handle



Telescopic Handle



Detachable Pillow



Safety Bandage



Wheel





## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901   
[www.chinasaikang.com](http://www.chinasaikang.com)

### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
Jiangsu Province,China.  
TEL: +86-512-58525850  
FAX: +86-512-58520713  
E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
17800 Castleton St Ste 180,  
City of Industry, Ca 91748  
Tel: +1 949648-2776

### SAIKANG (UK Branch)

Uk Healthward International Limited  
Owner(s) Address DEPT 114,  
196 High Road Wood Green,  
London, N22 8HH  
United Kingdom  
TEL: +44 2032399738

### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
Viraj Complex Unit 3 Mombasa Road,  
PO Box 14262-00800, Nairobi.  
TEL: +254 733-704 690

### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
X22-S03 England Cluster,  
International City, Dubai UAE  
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.