



SKB1C10

Electric
Stair Stretcher

SKW-B2



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098

DIMENSION PARAMETER

Full Size	1230*500*1600mm
Folded Size	1020*550*210mm

QUALITY CONTROL



STRUCTURAL DESIGN

Scientific design, ergonomics, simple structure, simple operation, reliable and secure. Widely used in the home for the old people, school security, hospitals, hotels etc for emergency evacuation.



FOLDABLE

Foldable stretcher, easy for carrying and storage



ELECTRIC CONTROL SYSTEM

Circuit is inside the pipe, waterproof, noiseless during operation. Compliance with EU safety standards. With LED Night light, LCD power display



WHEELS

With 4 wear-resistant rubber wheels, the chair can use as a wheelchair on the ground.

MECHANICAL FUNCTION

Can climb up and down stairs smoothly and steadily by rubber track.
Can move on the ground as manual wheelchair.
With Strong grip track, prevent glide or fall down when stop on the stairs.



POWDER COATING


Process Epoxy Painting, ASTM Testing Anti-Bacterial, Paint Thickness 0.12mm, Brightness 60°, Paint Can Resist 50kg Impact

0.12 | 60° | 50



RICH PARTS

LEADING BRAND OF CHINA MEDICAL FURNITURE



Folding Handle



Telescopic Handle



Detachable Pillow



Safety Bandage



Wheel





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.