



# SKE-120A

ELECTRIC  
DIALYSIS CHAIR



**SAIKANG** / Medical

Healthier  
your life

# SKE-120A

## ELECTRIC DIALYSIS CHAIR



### Technical parameters

▪ Size	1880×580mm
▪ Back	800×590mm
▪ Seat	490×590mm
▪ Leg	500×580mm
▪ Foot board	(230-360)×220mm
▪ Seat height	560mm
▪ Pillow	370×290mm
▪ Armrest	500×160mm
▪ Mattress thickness	60mm
▪ Back rest lifting	0° - 80°
▪ Leg rest lifting	20° - 90°

- Armrest swing out 0°-270°
- Voltage and frequency: 100-240 V · 50/60 Hz
- Plug as required (hospital grade)

### Technical configuration

<input checked="" type="checkbox"/> With CPR	
<input checked="" type="checkbox"/> Hand controller	1pc
<input checked="" type="checkbox"/> Power cord	1pc
<input checked="" type="checkbox"/> 3" silent caster	1set
<input checked="" type="checkbox"/> Small dining table	1pc



**Design**


European design, 4 sections ergonomics design.

**Steels**

Baosteel from Fortune Global 500.

**Plastic**

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.

**Welding**

Panasonic Robotic ensure 360° full smooth wedding.

**Electric Motor**

Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.

**Armrest**

Armrest updated, can be universal adjust. User can adjust any position he/she wants. Armrest is R120° concave curve, avoid donor's arm sliding, ensure the user be safety and comfortable.

**Leather**

Optional genuine leather or supebric cover,with artisanal making,comfortable hand feeling.

**Pillow**

Ergonomics design removable pillow.

**Table Board**

Adjustable ABS tray table, with antiskid convex edge, tank for cup and cellphone. USB/5V optional. Very convenient for donor.

**Dustproof Base Cover**

With dustproof base cover, protect motors and lines, easy to clean.

**IV Pole**

Full stainless steel IV pole, bearing 15 kg.



## ULTRAFIBER LEATHER COLOR TABLE



Off white 5509-1



Off white 5509-2



Beige 5509-3



Pink 5510-45




Light rose red 2892-15



Rose red 2201-30



Wine red 2209-22



Light blue 5510-43



Sky blue 5510-7



Dark lake blue 2201-40



Light grey 5509-14



Gray 5510-12



Lake blue 5509-27



Navy blue 2201-42



Light green 2210-88



Green 5510-20

Note: The above is the physical scan , there may be a certain color difference in kind , are to prevail in kind.



## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 📞 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue