



SKE001

ATTENDANT CHAIR



SKE001

FEATURES AND HIGHLIGHTS

Technical parameters


▪ Seat Size (LxWxH)	640×720×940 mm
▪ Bed Size (LxWxH)	640×1900×420 mm

Technical configuration

<input checked="" type="checkbox"/> Arm	1set
<input checked="" type="checkbox"/> 2" Castors	4pcs
<input checked="" type="checkbox"/> 3" Castors	2pcs
<input checked="" type="checkbox"/> Anti-slip Feet Cover	4pcs
<input checked="" type="checkbox"/> Pillow	1pc

Painting

1 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Headrest

The design of headrest is appropriate to protect cervical vertebra and make your sleep comfortable.



Design

Can pull out to be a bed, night chaperones are more comfortable.



Leather

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



Material

Baosteel from Fortune Global 500 Panasonic Robotic ensure 360° full smooth wedding.



Armrest

Frosted surface, anti-slip and durable.



Castors

High quality medical silent casters, 3" and 2" casters, anti-static electricity, protect the ground, waterproof and anti-bacteria.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue