



SKE004-1

TRANSFUSION CHAIR



SKE004-1

FEATURES AND HIGHLIGHTS

Technical parameters

- Seat Size (LxWxH) 670×800×1050 mm

Technical configuration

<input checked="" type="checkbox"/> Anti-slip Feet Cover	4pcs
<input checked="" type="checkbox"/> Wooden Armrest	2pcs
<input checked="" type="checkbox"/> IV Pole Previsions	2pcs
<input checked="" type="checkbox"/> IV Pole	1set



Design

Scientific design, ergonomics, simple structure, convenient operation, applied to infusion and patient rest.



Seat

The head has a sponge support, which can give the infusion person a sense of security and comfort during the infusion state.



Material

Plastic steel frame, strong stability, durable use, beautiful appearance, comfortable texture, convenient cleaning.



Armrest

The handrail is inlaid with painted wooden handrails. The appearance is beautiful, the atmosphere, and the effect of paint can better reflect the texture of the handrails and soft touch.



Back and Seat Plate

The back plate and seat plate are made of high-quality artificial leather and sponge, providing a safe, comfortable and warm infusion experience.



Foot Pad

Beautiful appearance, easy to use and has the functions of anti-slip, anti-collision, anti-shock, temperature resistance, noise reduction, etc.



IV Pole

Full stainless steel IV pole, bearing 15 kg.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue