



SKE091

BLOOD
DONATION CHAIR



SAIKANG / Medical

Healthier
your life

SKE091

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Seat Size (LxWxH)	560×1030×1200 mm
▪ Bed Size (LxWxH)	560×1650×470 mm
▪ Back-rest	0°-70°
▪ Leg Rest Down Fold	0°-90°

Technical configuration

<input checked="" type="checkbox"/> Arm	2pcs
<input checked="" type="checkbox"/> Anti-slip Feet Cover	4pcs

Foot Pad

Beautiful appearance, easy to use and has the functions of anti-slip, anti-collision, anti-shock, temperature resistance, noise reduction, etc.




Design

Through the button to achieve a key flat, seconds into a recliner, can lie down to sleep.



Chair Surface

Optional genuine leather or supebric cover, with artisanal making, comfortable hand feeling.



Material

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full smooth welding.



Armrest

Upgraded armrest, can be right-left, front-back adjustable. User can control universal position.



ULTRAFIBER LEATHER COLOR TABLE



Off white 5509-1



Off white 5509-2



Beige 5509-3



Pink 5510-45



Light rose red 2892-15



Rose red 2201-30



Wine red 2209-22



Light blue 5510-43



Sky blue 5510-7



Dark lake blue 2201-40



Light grey 5509-14



Gray 5510-12



Lake blue 5509-27



Navy blue 2201-42



Light green 2210-88



Green 5510-20

Note: The above is the physical scan , there may be a certain color difference in kind , are to prevail in kind.



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue