



SKE801

NURSE STOOL



SAIKANG / Medical

Healthier
your life

SKE801

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Stool surface	Φ350 mm
▪ Thickness of stool surface	80 mm
▪ Wheel height	490-630 mm
▪ Chassis	Φ400 mm



Design

This product adopts a quality steel pipe, corrosion resistance, no rust, no water, no moisture, safe and reliable.



Stool Surface

PU leather, high density sponge seat, soft and comfortable.



Gas Spring

Height can be adjusted, which is pneumatically controlled by handle underneath the seat.



Chassis

The chassis is made of high-quality stainless steel, it allows the feet to be placed more comfortably and sitting more comfortably.



Castors

Five nylon castors, convenient, wear-resistant, sliding silent.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue