



SKH046-14

OVERBED TABLE




SAIKANG / Medical

Healthier
your life

SKH046-14

FEATURES AND HIGHLIGHTS



Technical parameters

External size (LxWxH) 890×400×895-1195mm

Technical configuration

- | | |
|--|-----|
| <input checked="" type="checkbox"/> Air Spring for Height Adjustment | 1pc |
| <input checked="" type="checkbox"/> Phenolic Resin Table Top | 1pc |

Material

U type base is made of powder coated steel, can be used in combination with bedside cabinets.



Table Top

Table top is made of HDPE, easy clean, high quality environmental material, light weight, corrosion resistance, impact resistance, easy to clean and disinfect, durable for use. With groove to put cups and pen.



Height Adjustment

Height adjustment by gas spring, used single-handed according to the convenience of the user.



Casters

Easy to install, more stable and firmer. High quality medical caster, 1.5 inch / 2 inch mute wheel.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue