



# SKR-AT625-1

## ANESTHESIA TROLLEY

Registration code: SKM-B



**SAIKANG** / Medical

Healthier  
your life

# SKR-AT625-1

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 625×475×930 mm

### Technical configuration

<input checked="" type="checkbox"/> Sliding Side Shelf	1pc
<input checked="" type="checkbox"/> Dust Baskets	1pc
<input checked="" type="checkbox"/> Utility Container	1pc
<input checked="" type="checkbox"/> Central Lock	1set
<input checked="" type="checkbox"/> Trash Can	2pcs
<input checked="" type="checkbox"/> Anesthesia Frame	1pc



#### Powder Coated Column

Make the structure of the trolley more stable and firm.




#### Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley.



#### Drawers and Cabinet

Has two drawers, different size. With doors, can put more things.



#### Sliding Tray

Laterally extendable worktop.



#### File Box

Side attached file box, for storage of paper work.



#### Debris Basket

It is easy to put some files about patient records.



#### Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



#### Trash Can

There are two ABS dirt buckets of the car body, different colors can distinguish different kinds of dirt.




#### Drawer Card Holder

card insert from downside, easy keep clean.



#### Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.



#### Anesthesia Frame

Anesthesia frame inserted into main body column.



#### Heavy Duty Construction

Wheels assembled on the same vertical line as trolley four columns, more stable.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue