



SKR-ET625

EMERGENCY TROLLEY

Registration code: SKM-B



SAIKANG / Medical

Healthier
your life

SKR-ET625

FEATURES AND HIGHLIGHTS

Technical parameters

- Size (LxWxH) 625×475×930 mm

Technical configuration

<input checked="" type="checkbox"/> IV pole	1 pc.
<input checked="" type="checkbox"/> Sliding side shelf	1 pc.
<input checked="" type="checkbox"/> CPR board	1 pc.
<input checked="" type="checkbox"/> Dust basket	1 pc.
<input checked="" type="checkbox"/> Utility container	1 pc.
<input checked="" type="checkbox"/> Trash can	2 pcs.
<input checked="" type="checkbox"/> Castors	4 pcs.
<input checked="" type="checkbox"/> Central lock	1 set
<input checked="" type="checkbox"/> Power outlet	1 pc.
<input checked="" type="checkbox"/> Defibrillator shelf	1 pc.
<input checked="" type="checkbox"/> Oxygen bottle holder	1 pc.



Powder Coated Column

Make the structure of the trolley more stable and firm.



Table Top

Stainless steel side rail, avoid possible fall down during movement of trolley.



Drawers

1-3 layer small drawer, fourth floor middle drawer, fifth floor large drawer, can put different sizes of medicine boxes and other medical materials, and equipped with central control lock.



Sliding Tray

Laterally extendable worktop.



File Box

Side attached file box, for storage of paper work.



Debris Basket

It is easy to put some files about patient records.



Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.

**Trash Can**

There are two ABS dirt buckets of the car body, different colors can distinguish different kinds of dirt.

**Drawer Card Holder**

card insert from downside, easy keep clean.

**Castor**

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.

**Heavy Duty Construction**

Wheels assembled on the same vertical line as trolley four columns, more stable.

**Power Cord**

Power cord & sockets detachable.

**CPR Board**

Used to save apatient in an emergency.

**Oxygen Cylinder Holder**

Inserted oxygen bottle holder, belts fixed on trolley main body.

**Defibrillator Holder**

The plate surface is made of virgin ABS material. The position of the IV pole and the defibrillator can be exchanged for easy use.

**IV Pole**

Hidden height adjustable IV pole, 4 hooks, dia 25mm.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue