

SKS008

ABS+STEEL
BEDSIDE TABLE



SKS008

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 500x510x705 mm



Design

The bedside table is designed with complete functions, beautiful appearance, not easy to deform, etc., and it is easy to clean and take care of, and it is durable.



Material

The cabinet body is made of high-quality steel sprayed plastic. The countertops, drawers, drawers and cabinet doors are all made of ABS, 100% pure new material injection molding, non-toxic and harmless, sturdy, rust and corrosion resistant, easy to clean, and have a long service life.



Towel Hook

There are concealed towel hooks on both sides of the cabinet, which can be pulled out when in use and hidden when not in use, saving space, strong, rust-resistant, and easy to clean.



Cabinet Body

There is a drawer as a whole, and there is a pull-out dining board between the drawer and the table top; a cabinet with two layers inside and large storage space.



Castors

Silent castors, easy to install, strong and stable, easy to transfer, anti-winding.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue