



# SKB041-1

## PATIENT MANUAL STRETCHER

Registration code: SKM-F





**SAIKANG** / Medical

Healthier  
your life

# SKB041-1

## FEATURES AND HIGHLIGHTS



### Technical parameters

▪ External size (LxW):	1930x640 mm
▪ Table height adjustment:	540-840 mm
▪ Back plate adjustment:	0-75° ( $\pm 10^\circ$ )

### Technical configuration

<input checked="" type="checkbox"/> Crank system	1 set
<input checked="" type="checkbox"/> Back plate air spring	1 set
<input checked="" type="checkbox"/> ABS side guardrail	1 set
<input checked="" type="checkbox"/> 6" castors with central locking system	4 pcs.
<input checked="" type="checkbox"/> Directional wheel	1 pc.
<input checked="" type="checkbox"/> Central locking pedal	2 pcs.
<input checked="" type="checkbox"/> Directional pedal	2 pcs.
<input checked="" type="checkbox"/> IV pole	1 set
<input checked="" type="checkbox"/> IV pole prevision	2 pcs.
<input checked="" type="checkbox"/> Oxygen bottle holder	1 pc.
<input checked="" type="checkbox"/> Transfer mattress	1 pc.



**Metal material**

Baosteel from Fortune Global 500  
Panasonic Robotic ensure 360° full  
smooth welding.

**Bearing**

Robust tube in weight bearing  
structure.

**Painting**

11 process epoxy painting, ASTM  
testing anti-bacterial, paint thickness  
0.12 mm, brightness 60°. Painting can  
resist 50 kg impact.

**Protection**

Cover all metal screws with rubber  
cap, protect screws from possible  
corrosion.  
Metal dust protection for crank  
system, more durable than plastic  
cover.

**Plastic**

Environmental plastic can be  
sterilized in 100°C, tensile upto 30  
MPa.

**Back lift**

Back lift can be up to 80°, making  
patient feel more comfortable when  
needed.

**Double side castors**

Four double side castors of 6 inch  
diameter, simultaneously lockable by  
the step pedal at both sides.  
TPR tire no worn out after running  
30 km, save anti-winding hard shell,  
united forming without bolts.

**Directional center wheel**


Can be engaged to facilitate direction  
control, enabling users to navigate  
the stretcher easily in hectic and  
dynamic environments.

**Telescopic IV pole**

Full stainless steel IV pole, bearing  
15 kg.

**Oxygen bottle holder**

Adjustable vertical oxygen bottle  
holder. Load capacity: 15 kg.

**Stainless steel handle**

Stainless steel handle, more durable.

**Matress**

Water proof matress with tow belts.







## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi,  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai,  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue