

# SKH038-2

## MAYO TROLLEY

Registration code: SKM-B



# SKH038-2

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 500×350×(870-1300) mm
- Tray size (LxWxH) 350×500×48 mm

### Technical configuration

- Castor 4 pcs
- Tray 1 pc



#### Material

Φ25mm, Φ16mm stainless steel round pipe, the whole frame used high quality S.S. with the Mirror surface polishing.



#### Craftsmanship

Imported intelligent robot is used for welding forming, high precision, no blind spots, smooth surface, no burrs, beautiful and strong, high bearing capacity.



#### Tray

Size of stainless steel tray:  
L350\*W500\*H48mm.



#### Castor

3 inch covered castors (diagonal brake), flexible and light to drive, smooth steering, wear-resistant, two with brake, can be easily braked in any state.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue