



SKB041-3S

HYDRAULIC PATIENT TROLLEY

Registration code: SKM-G





SAIKANG / Medical

Healthier
your life

SKB041-3S

FEATURES AND HIGHLIGHTS



Technical parameters

- External size (LxW): 1950x650 mm
- Table height adjustment: 620-950 mm
- Back plate adjustment: 0-55° ($\pm 10^\circ$)


Technical configuration

- Hydraulic Pump 1 pc
- Air Spring 1 pc
- Stainless steel guardrail 1 set
- Push handle 1 set
- Wall bumper wheel 4 pcs
- Foot pedal 2 pcs
- 5" Castors with brake 4 pcs
- Oxygen cylinder rack 1 pc
- Debris basin 1 pc
- IV pole 1 pc
- Mattress with thickness of 60mm 1 pc



Metal material

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full
smooth welding.

**Painting**

11 process epoxy painting, ASTM
testing anti-bacterial, paint thickness
0.12 mm, brightness 60°. Painting can
resist 50 kg impact.

**Multifunctional patient cart**

Multi-functional trolley, which can
be separated into two parts: stretcher
and trolley.
The stretcher with four wheels can be
used and moved individually.

**Bumper wheel**

Bumper wheel detachable, max.
protect the bed and patient during
transportation.

**Guardrail**

Foldable stainless steel guardrail, no
rust, easy to clean.

**Hydraulic pumps**

American brand hydraulic pump,
providing strong support for the
functions of height adjustment.

**Casters**


Four covered casters of 5 inch
diameter, all with brakes.

**IV pole**

Foldable telescopic stainless steel IV
pole, bearing 15kg.

**Oxygen bottle holder**

Oxygen bottle holder, load capacity:
15KG.

**Stainless steel basin**

Equipped with stainless steel sewage
basin and other various accessories.

**Mattress**

Waterproof mattress, convenient for
transfer patient.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue