



R000w

ONE FUNCTION MANUAL BED

Registration code: SK-CD1




SAIKANG Medical

Healthier
your life

R000w

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxWxH)	2000x1000x490 mm
▪ Mattress platform	1950x950 mm
▪ Safe working load	200 kg
▪ Back-rest adjustment	0-75° (±10°)



Technical configuration

<input checked="" type="checkbox"/> Powder coated Bed Ends	1set
<input checked="" type="checkbox"/> IV Pole	1pc
<input checked="" type="checkbox"/> IV Pole Prevision	1p

Optional

<input type="checkbox"/> 5" Noiseless Castors	
<input type="checkbox"/> Integrated detachable guardrail	



Steels | Painting


Baosteel from Fortune Global500.
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12 mm, brightness 60°, paint can resist 50kg impact.

**Welding**

Panasonic robotic ensure, 360° full smooth welding

**Bed Surface**

Reticulated bed platform.

**Head & Foot Board**

Powder coated head&foot board, simple installation, Removable.

**Handle**


Equipped with a handle to manually control the angle of back lift.

**IV Pole**

Full stainless steel IV pole, bearing 15 kg.

**Optional**

Integrated detachable guardrail.

**Optional**

5" noiseless castors.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue