



SKB039(E)

AMBULANCE STRETCHER TROLLEY

Registration code: SKC-3F5



SAIKANG / Medical

Healthier
your life

SKB039(E)

FEATURES AND HIGHLIGHTS

Technical parameters

▪ Full size (LxWxH)	1900x540x800mm
▪ Lowest size (LxWxH)	1900x540x200mm
▪ Back-rest lifting	0-85°
▪ Net weight	45kg
▪ Load bearing	≤159kg

Side rail

Side rail is foldable.



Safety lock

With safety lock, it must be locked up when transport the patients, and opened out when get on the ambulance.



Caster

5' silence caster. Advanced protection universal brake wheel, easy to move, simple to operate, reliable performance.



Function

Can be folded into a chair, the angle of the chair can be adjusted in five gear.



Fixing devices

Can be locked by the fixing devices when on the ambulance car.



Structure

Premium alluminium alloy structure, more stable, but light weight.



Leather / Safety belts

Artificial leather, bright appearance, soft texture, abrasion resistance, folding resistance, acid and alkali resistance. With 2 safety belts, make patient safer.



Sponge

Flame retardant, 5 cm thickness, high-density foam, density 30 kg/m³.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 📞 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue